

Supplementary Information

Table S3

Information on the primer sequences for qRT-PCR analysis in this study

Gene	Forward (5'-3')	Reverse (5'-3')
miR-9	CGGTCTTGTTATCTAGCTGTATG	RR716, SYBR® PrimeScript™ miRNA RT-PCR Kit, TaKaRa
miR-9-1	CAGAGAAGGGCAGTGGAGAC	ACGACAGAGACCGAAAAAGG
miR-9-2	CTAACGCTGCCGGAGATTAC	TACTTGC CGCGCTTAAGATT
miR-9-3	GCGCTCGAGGCTCTCTAAG	GAGGGGATGGACAGACACAC
COL18A1	GCCATACTTCCTGTATAGTTC	AGGTGAATGTCAGCATCA
THBS2	GATTAGAACAAACAGGCATAG	GTAGGACAAGAGGGAGAGTTC
PTCH1	AAGTTGGAGGACGAGTAA	GACATTAGCACCTTCTTCTT
PHD3	TTCTAACGACTGCTCTATT	CCTGTCTTCATTGTGCTTAC
HIF-1α	CCGAGGAAGAACTATGAA	GTTGGTTACTGTTGGTATC
VEGF	AGCCTTGCCTTGCTGCTCTA	GTGCTGGCCTTGGTGAGG
MYC	CTCCTGGCAAAAGGT CAGAG	TCGGTTGTTGCTGATCTGTC
OCT4	GTATT CAGCCAAACGACCATCT	GCTTCCTCCACCCACTTCT
U6	CTCGCTTCGGCAGCACA	AACGCTTCACGAATTGCGT
β-actin	TGGCATCCACGAAACTACC	GTGTTGGCGTACAGGTCTT