



CHEMOTECHNIQUE DIAGNOSTICS

Compound	Conc. Veh. %(w/w)	Art. No.
European Baseline Series		S-1000
1. Potassium dichromate	0.5 pet	P-014A
2. p-PHENYLENEDIAMINE (PPD)	1.0 pet	P-006
3. Thiuram mix	1.0 pet	Mx-01
-Dipentamethylenethiuram disulfide	0.25	D-019
-Tetraethylthiuram disulfide (TETD)	0.25	T-002
-Tetramethylthiuram disulfide (TMTD)	0.25	T-005
-Tetramethylthiuram monosulfide (TMTM)	0.25	T-006
4. Neomycin sulfate	20.0 pet	N-001
5. Cobalt(II)chloride hexahydrate	1.0 pet	C-017A
6. Benzocaine	5.0 pet	B-004
7. Nickel(II)sulfate hexahydrate	5.0 pet	N-002A
8. Clioquinol	5.0 pet	C-015
9. COLOPHONIUM	20.0 pet	C-020
10. Paraben mix	16.0 pet	Mx-03C
-BUTYLPARABEN	4.0	B-020
-ETHYLPARABEN	4.0	E-010
-METHYLPARABEN	4.0	M-012
-PROPYLPARABEN	4.0	P-020
11. N-Isopropyl-N-phenyl-4-phenylenediamine (IPPD)	0.1 pet	I-004
12. LANOLIN ALCOHOL	30.0 pet	W-001
13. Mercapto mix	2.0 pet	Mx-05A
-N-Cyclohexyl-2-benzothiazolesulfenamide	0.5	C-023
-2-Mercaptobenzothiazole (MBT)	0.5	M-003
-Dibenzothiazyl disulfide (MBTS)	0.5	D-003
-2-(4-Morpholinylmercapto)benzothiazol (MOR)	0.5	M-016
14. Epoxy resin, Bisphenol A	1.0 pet	E-002
15. MYROXYLON PEREIRAE RESIN*	25.0 pet	B-001
16. 4-tert-Butylphenolformaldehyde resin (PTBP)	1.0 pet	B-024
17. 2-Mercaptobenzothiazole (MBT)	2.0 pet	M-003A
18. FORMALDEHYDE	2.0 aq	F-002B
19. Fragrance mix I*	8.0 pet	Mx-07
-AMYL CINNAMAL	1.0	A-014

* Emulsifier: SORBITAN SESQUIOLEATE 5%

...world leader in patch testing

Compound		Conc. Veh. %(w/w)	Art. No.
	-CINNAMYL ALCOHOL	1.0	C-013
	-CINNAMAL	1.0	C-014
	-EUGENOL	1.0	E-016
	-GERANIOL	1.0	G-001
	-HYDROXYCITRONELLAL	1.0	H-008
	-ISOEUGENOL	1.0	I-002
	-Oakmoss absolute	1.0	O-001
20.	Sesquiterpene lactone mix	0.1 pet	Mx-18
	-Alantolactone	0.033	A-003
	-Costunolide	0.033	C-039
	-Dehydrocostus lactone	0.033	D-056
21.	QUATERNIUM-15	1.0 pet	C-007A
22.	2-Methoxy-6-n-pentyl-4-benzoquinone	0.01 pet	M-008
23.	METHYLISOTHIAZOLINONE + METHYLCHLOROISOTHIAZOLINONE	0.01 aq	C-009A
24.	Budesonide	0.01 pet	B-033B
25.	Tixocortol-21-pivalate	0.1 pet	T-031B
26.	METHYLDIBROMO GLUTARONITRILE	0.5 pet	D-049E
27.	Fragrance mix II	14.0 pet	Mx-25
	-Hexyl cinnamic aldehyde	5.0	H-025
	-COUMARIN	2.5	C-038
	-FARNESOL	2.5	F-004
	-Lyrall	2.5	L-003
	-CITRAL	1.0	C-036
	-CITRONELLOL	0.5	C-037
28.	Lyrall	5.0 pet	L-003
29.	METHYLISOTHIAZOLINONE	0.2 aq	M-035B
Revised January 2014			

It is strongly recommended to make an additional reading of the test on day 7.

References:

1. M. Isaksson, F.M. Brandao, M.Bruze, A.Goossens. Recommendation to include budesonide and tixocortol privalate in the european standard series. *Contact dermatitis* 2000, 43, 41-42.
2. Maureen J. Jonker & Derk P. Bruynzeel. The outcome of an additional test reading on days 6 or 7. *Contact dermatitis* 2000, 42, 330-335.
3. Recommendation to include Fragrance mix II and hydroxyisohexyl 3-cyclohexene car-



CHEMOTECHNIQUE DIAGNOSTICS

Compound

Conc. Veh.
%(w/w)

Art. No.

boxaldehyde (Lyrat®) in Contact Dermatitis Volume 58, Issue 3, March 2008, Pages: 129–133, Magnus Bruze, Klaus Ejner Andersen and An Goossens, on behalf of the ESCD and EECDRG

Various national baseline series developed by the national contact dermatitis research groups are also available on request.

International Comprehensive Baseline Series

ICB-1000

1.	Benzocaine*	5.0	pet	B-004
2.	2-Mercaptobenzothiazole (MBT)	1.0	pet	M-003B
3.	COLOPHONIUM*	20.0	pet	C-020
4.	p-PHENYLENEDIAMINE (PPD)*	1.0	pet	P-006
5.	IMIDAZOLIDINYL UREA	2.0	pet	I-001A
6.	CINNAMAL	1.0	pet	C-014
7.	Amerchol L 101	50.0	pet	A-004
8.	Carba mix	3.0	pet	Mx-06
	-1,3-Diphenylguanidine	1.0		D-022
	-ZINC DIBUTYLDITHIOCARBAMATE (ZBC)	1.0		Z-002
	-Zinc diethyldithiocarbamate (ZDC)	1.0		Z-003
9.	Neomycin sulfate*	20.0	pet	N-001
10.	Thiuram mix*	1.0	pet	Mx-01
	-Dipentamethylenethiuram disulfide	0.25		D-019
	-Tetraethylthiuram disulfide (TETD)	0.25		T-002
	-Tetramethylthiuram disulfide (TMTD)	0.25		T-005
	-Tetramethylthiuram monosulfide (TMTM)	0.25		T-006
11.	Clobetasol-17-propionate	1.0	pet	C-028
12.	Ethylenediamine dihydrochloride	1.0	pet	E-005
13.	Epoxy resin, Bisphenol A*	1.0	pet	E-002
14.	QUATERNIUM-15*	2.0	pet	C-007B
15.	4-tert-Butylphenolformaldehyde resin (PTBP)*	1.0	pet	B-024
16.	Mercapto mix	1.0	pet	Mx-05B
	-N-Cyclohexyl-2-benzothiazolesulfenamide	0.25		C-023
	-2-Mercaptobenzothiazole (MBT)	0.25		M-003

* Also present in European Baseline Series

...for the diagnosis of contact allergy

Compound	Conc. Veh. %(w/w)	Art. No.
-Dibenzothiazyl disulfide (MBTS)	0.25	D-003
-2-(4-Morpholinylmercapto)benzothiazol (MOR)	0.25	M-016
17. N-Isopropyl-N-phenyl-4-phenylenediamine (IPPD)*	0.1 pet	I-004
18. Potassium dichromate	0.25 pet	P-014B
19. MYROXYLON PEREIRAE RESIN ^{*,**}	25.0 pet	B-001
20. Nickel(II)sulfate hexahydrate	2.5 pet	N-002B
21. DIAZOLIDINYL UREA	1.0 pet	D-044C
22. TOCOPHEROL	100	T-036
23. Bacitracin	20.0 pet	B-032B
24. Mixed dialkyl thiourea	1.0 pet	Mx-24
-N,N'-Dibutylthiourea	0.5	D-038
-N,N'-Diethylthiourea	0.5	D-039
25. DISPERSE ORANGE 3	1.0 pet	D-032
26. Paraben mix	12.0 pet	Mx-03A
-BUTYLPARABEN	3.0	B-020
-ETHYLPARABEN	3.0	E-010
-METHYLPARABEN	3.0	M-012
-PROPYLPARABEN	3.0	P-020
27. METHYLDIBROMO GLUTARONITRILE*	0.5 pet	D-049E
28. Fragrance mix I ^{*,**}	8.0 pet	Mx-07
-AMYL CINNAMAL	1.0	A-014
-CINNAMYL ALCOHOL	1.0	C-013
-CINNAMAL	1.0	C-014
-EUGENOL	1.0	E-016
-GERANIOL	1.0	G-001
-HYDROXYCITRONELLAL	1.0	H-008
-ISOEUGENOL	1.0	I-002
-Oakmoss absolute	1.0	O-001
29. GLUTARAL ^{**}	0.5 pet	G-003B
30. 2-BROMO-2-NITROPROPANE-1,3-DIOL	0.5 pet	B-015B
31. Sesquiterpene lactone mix [*]	0.1 pet	Mx-18
-Alantolactone	0.033	A-003
-Costunolide	0.033	C-039
-Dehydrocostus lactone	0.033	D-056

* Also present in European Baseline Series

** Emulsifier: SORBITAN SESQUIOLEATE 5%



CHEMOTECHNIQUE DIAGNOSTICS

Compound	Conc. Veh. %(w/w)	Art. No.
32. THIMEROSAL	0.1 pet	T-007
33. PROPOLIS	10.0 pet	P-022
34. BENZOPHENONE-3	10.0 pet	H-014C
35. CHLOROXYLENOL (PCMX)	1.0 pet	C-010B
36. Ethyleneurea, melamine formaldehyde mix**	5.0 pet	Mx-16
-Dimethylol dihydroxy ethylene urea	4.0	D-012
-Melamine formaldehyde	1.0	M-001
37. 2-tert-Butyl-4-methoxyphenol (BHA)	2.0 pet	B-022
38. Gold(I)sodium thiosulfate dihydrate	0.5 pet	G-005A
39. Ethyl acrylate	0.1 pet	E-004
40. GLYCERYL THIOGLYCOLATE	1.0 pet	G-004
41. Toluenesulfonamide formaldehyde resin	10.0 pet	T-010
42. Methyl methacrylate	2.0 pet	M-013
43. Cobalt(II)chloride hexahydrate*	1.0 pet	C-017A
44. Tixocortol-21-pivalate*	0.1 pet	T-031B
45. Budesonide*	0.01 pet	B-033B
46. COCAMIDE DEA	0.5 pet	C-019
47. TRIETHANOLAMINE	2.0 pet	T-016
48. Hydrocortisone-17-butyrate	1.0 pet	H-021B
49. Tea Tree Oil oxidized	5.0 pet	T-035B
50. Fragrance mix II*	14.0 pet	Mx-25
-Hexyl cinnamic aldehyde	5.0	H-025
-COUMARIN	2.5	C-038
-FARNESOL	2.5	F-004
-Lyr al	2.5	L-003
-CITRAL	1.0	C-036
-CITRONELLOL	0.5	C-037
51. Disperse Yellow 3	1.0 pet	D-036
52. BENZYL SALICYLATE	10.0 pet	B-010B
53. DECYL GLUCOSIDE***	5.0 pet	D-065
54. METHYLISOTHIAZOLINONE*	0.2 aq	M-035B
55. 2-Hydroxyethyl methacrylate	2.0 pet	H-010
56. DMDM HYDANTOIN†***	1.0 pet	D-047B

* Also present in European Baseline Series

** Emulsifier: SORBITAN SESQUIOLEATE 5%

*** Emulsifier: SORBITAN SESQUIOLEATE 1%

...the trusted name in patch testing

Compound	Conc. Veh. %(w/w)	Art. No.
57. CANANGA ODORATA OIL	2.0 pet	Y-001
58. BENZYL ALCOHOL	10.0 sof	B-008B
59. ISOPROPYL MYRISTATE	20.0 pet	I-003
60. TRICLOSAN	2.0 pet	T-014
61. Desoximetasone	1.0 pet	D-057
62. POLYSORBATE 80	5.0 pet	P-013
63. IODOPROPYNYL BUTYLCARBAMATE	0.2 pet	I-008C
64. 2-n-Octyl-4-isothiazolin-3-one	0.1 pet	O-004
65. Disperse Blue mix 106/124	1.0 pet	Mx-26
-Disperse Blue 106	0.5	D-040
-Disperse Blue 124	0.5	D-041
66. Compositae mix II	5.0 pet	Mx-29A
-ANTHEMIS NOBILIS EXTRACT	1.2	C-029
-CHAMOMILLA RECUTITA EXTRACT	1.2	C-051
-ACHILLEA MILLEFOLIUM EXTRACT	1.0	A-025
-TANACETUM VULGARE EXTRACT	1.0	T-033
-ARNICA MONTANA EXTRACT	0.5	A-024
-Parthenolide	0.1	P-029
67. Lidocaine	15.0 pet	L-002B
68. Fusidic acid sodium salt	2.0 pet	F-003
69. Dibucaine hydrochloride	2.5 pet	D-005B
70. Benzoylperoxide	1.0 pet	B-007
71. ISOAMYL p-METHOXYCINNAMATE	10.0 pet	I-009
72. Lyr ^{al} *	5.0 pet	L-003
73. ETHYLHEXYL SALICYLATE	5.0 pet	O-007A
74. BENZALKONIUM CHLORIDE	0.1 aq	B-027
75. Amidoamine	0.1 aq	A-029
76. COCAMIDOPROPYL BETAINE	1.0 aq	C-018
77. FORMALDEHYDE*	2.0 aq	F-002B
78. METHYLISOTHIAZOLINONE + METHYLCHLOROISOTHIAZOLINONE*	0.01 aq	C-009A
79. PROPYLENE GLYCOL	30.0 aq	P-019B
80. Dimethylol dihydroxy ethylene urea	4.5 aq	D-012
81. Hydroperoxides of Linalool	1.0 pet	H-031
82. Hydroperoxides of Limonene	0.3 pet	H-032

* Also present in European Baseline Series

Revised January 2014

1.	N-Isopropyl-N-phenyl-4-phenylenediamine (IPPD)*	0.1	pet	I-004
2.	GLUTARAL**	0.2	pet	G-003A
3.	DISPERSE ORANGE 3	1.0	pet	D-032
4.	Acid yellow 36	1.0	pet	A-019
5.	Hydroquinone monobenzylether	1.0	pet	H-019
6.	Thiuram mix*	1.0	pet	Mx-01
	-Dipentamethylene	0.25		D-019
	-Tetraethylthiuram	0.25		T-002
	-Tetramethylthiuram	0.25		T-005
	-Tetramethylthiuram monosulfide (TMTM)	0.25		T-006
7.	Potassium dichromate*	0.5	pet	P-014A
8.	4-tert-Butylphenolformaldehyde resin (PTBP)*	1.0	pet	B-024
9.	p-PHENYLENEDIAMINE (PPD)*	1.0	pet	P-006
10.	Nickel(II)sulfate hexahydrate*	5.0	pet	N-002A
11.	COLOPHONIUM*	20.0	pet	C-020
12.	FORMALDEHYDE*	2.0	aq	F-002B
13.	N,N'-Diphenylthiourea (DPTU)	1.0	pet	D-025
14.	2-Mercaptobenzothiazole (MBT)*	2.0	pet	M-003A
15.	N,N'-Diethylthiourea	1.0	pet	D-039
16.	1,3-Diphenylguanidine	1.0	pet	D-022
17.	N,N'-Dibutylthiourea	1.0	pet	D-038
18.	Epoxy resin, Bisphenol A*	1.0	pet	E-002
19.	Dodecyl mercaptan	0.1	pet	D-043
20.	METHYLISOTHIAZOLINONE + METHYLCHLORO-ISOTHIAZOLINONE	0.02	aq	C-009B
21.	4-Aminoazobenzene	0.25	pet	A-005
22.	2-n-Octyl-4-isothiazolin-3-one	0.1	pet	O-004
23.	4,4'-Dithiodimorpholine	1.0	pet	D-054

Revised January 2014

* Also present in European Baseline Series

** Emulsifier: SORBITAN SESQUIOLEATE 5%