

Table S1. Overview of tet-resistance genes surveyed in the guts of *G. mellonella* larvae.

gene name	gene product	reference
<i>tetA</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetB</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetC</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetD</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetE</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetG</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetH</i>	efflux pump	Fan <i>et al.</i> , 2007 [1]
<i>tetJ</i>	efflux pump	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetK</i>	efflux pump	Chopra and Roberts, 2001 [4]
<i>tetL</i>	efflux pump	Tian <i>et al.</i> , 2012 [5]
<i>tetY</i>	efflux pump	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetZ</i>	efflux pump	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tet30</i>	efflux pump	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetBP</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetM</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetO</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetQ</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetS</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetT</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tetW</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tet32</i>	ribosomal protection protein	Aminov <i>et al.</i> , 2002, 2004 [2, 3]
<i>tet34</i>	efflux pump	Szczepanowski <i>et al.</i> , 2009 [6]
<i>tetX</i>	tetracycline destructase	Bartha <i>et al.</i> , 2011 [7]
<i>tetX2</i>	tetracycline destructase	Bartha <i>et al.</i> , 2011 [7]

References:

1. Fan W, Hamilton T, Webster-Sesay S, Nikolich MP, Lindler LE: **Multiplex real-time SYBR Green I PCR assay for detection of tetracycline efflux genes of Gram-negative bacteria.** *Mol Cell Probes* 2007, **21**(4):245-256.
2. Aminov RI, Chee-Sanford JC, Garrigues N, Mehboob A, Mackie RI: **Detection of tetracycline resistance genes by PCR methods.** *Methods Mol Biol* 2004, **268**:3-13.
3. Aminov RI, Chee-Sanford JC, Garrigues N, Teferedegne B, Krapac IJ, White BA, Mackie RI: **Development, validation, and application of PCR primers for detection of tetracycline efflux genes of gram-negative bacteria.** *Appl Environ Microbiol* 2002, **68**(4):1786-1793.
4. Chopra I, Roberts M: **Tetracycline antibiotics: mode of action, applications, molecular biology, and epidemiology of bacterial resistance.** *Microbiol Mol Biol Rev* 2001, **65**(2):232-260.
5. Tian B, Fadhil NH, Powell JE, Kwong WK, Moran NA: **Long-term exposure to antibiotics has caused accumulation of resistance determinants in the gut microbiota of honeybees.** *mBio* 2012, **3**(6).
6. Szczepanowski R, Linke B, Krahn I, Gartemann KH, Gutzkow T, Eichler W, Puhler A, Schluter A: **Detection of 140 clinically relevant antibiotic-resistance genes in the plasmid metagenome of wastewater treatment plant bacteria showing reduced susceptibility to selected antibiotics.** *Microbiology* 2009, **155**(Pt 7):2306-2319.
7. Bartha NA, Soki J, Urban E, Nagy E: **Investigation of the prevalence of tetQ, tetX and tetX1 genes in *Bacteroides* strains with elevated tigecycline minimum inhibitory concentrations.** *International journal of antimicrobial agents* 2011, **38**(6):522-525.